


CollectEject™ Swabs



Features

- DNA yield up to 5 µg from oral sample
- Collect an oral sample in 1 minute
- Sterile swab with multi-layer, ejectable swab head for maximum sample absorption and extraction
- Samples are stabilized at ambient temperature

CollectEject	Catalog #	Yield	Swab Head Length	Sample Type
Oral	G-2092-000	Up to 5 µg DNA	16mm	Oral Samples
Large Oral	G-2007-000	Up to 10 µg DNA	22mm	Oral Samples
Touch DNA	G-2070-100	N/A	16mm	DNA from surfaces
Blade DNA	G-2077-100	N/A	21.5mm	DNA from surfaces



Publication

PiBarreck M, Moreira F, Foley M. Evaluation of Maximum DNA Yield from a New Noninvasive Buccal Collection Device Following Various Extraction Protocols. *Genet Test Mol Biomarkers*. 2021 May;25(5):376-380. doi: 10.1089/gtmb.2020.0341. Epub 2021 Apr 29. PMID: 33926219

